



Fann Instrument Company

Included Parts

100071979
SLURRY CUP ASSY - AUTOMATED HPHT

Part Number	Description	Qty
100004335	PLUG BASE SLURRY CUP CONSISTOMETER	1
100010310	GASKET CUP PLUG	1
100007921	BASE SLURRY CUP f/AUTO HPHT CONSISTOMTER	1
100010309	GASKET CUP BASE HPHT CONSISTOMETER	1
100007920	SLEEVE SLURRY CUP	1
100004333	SHAFT SLURRY CUP	1
100004313	PADDLE SLURRY CUP	1
100021371	PIN ROLL 1/16"DIA. X 1/2" LONG STAINLESS	1
100007922	COLLAR DIAPHRAGM SLURRY CUP	1
100013735	DIAPHRAGM NEOPRENE SLURRY CUP	1
100007917	SUPPORT DIAPHRAGM SLURRY CUP	1
100007916	RING LOCK .625 LONG SLURRY CUP	1
100007918	DISC DRIVE SLURRY CUP	1
100007919	BAR DRIVE SLURRY CUP HPHT CONSISTOMETER	1
100072504	PIN, SHEAR, SLURRY CUP, HP-HT CONSISTOMETER	1
100002475	PIN ROLL 1/4 X 1/2 STAINLESS	2
207764	10-32 X 3/8 HSSS BOPL	1