



Fann Instrument Company

Included Parts

100072538

SLURRY CUP ASSY THIXOTROPIC CEMENT HPHT CONSISTOMETER

Part Number	Description	Qty
100072543	BAFFLE RETAINING THIXOTROPIC CONSISTOMETER	1
100072539	SLEEVE SLURRY CUP THIXO CONSISTOMETER	1
100072542	SLEEVE ASSEMBLY SLURRY CUP	1
100072541	RING RETAINING SIDE BLADES	1
100072540	BASE SLURRY CUP f/THIXOTROPIC CONSISTOMETER	1
100004335	PLUG BASE SLURRY CUP CONSISTOMETER	1
100072546	DIAPHRAGM BUNA-N SLURRY CUP	1
100072536	RING, LOCK, SLURRY, CUP	1
100072547	PADDLE, SLURRY CUP THIXO CONSISTOMETER	1
100072548	BLADE ASSEMBLY SIDE THIXO CONSISTOMETER	1
100021431	PIN, ROLL, 1/16 X 3/4, STAINLESS STEEL	1
100007922	COLLAR DIAPHRAGM SLURRY CUP	1
100072544	RING, SPACER, TOP BAFFLE, SLURRY CUP	1
100072545	BAFFLE TOP SLURRY CUP THIXO CONSISTOMETER	1
100004333	SHAFT SLURRY CUP	1
100004336	HUB SLURRY CUP THIXOTROPIC CONSISTOMETER	1
100004337	CAP HUB SLURRY CUP THIXO CONSISTOMETER	1
100004338	RING PACKING HUB SLURRY THIXO CONSISTOMETER	1
100002427	ORING 3/8 X 1/4 X 1/16 TEFLON	2
100010309	GASKET CUP BASE HPHT CONSISTOMETER	1
100007918	DISC DRIVE SLURRY CUP	1
100007919	BAR DRIVE SLURRY CUP, HPHT CONSISTOMETER	1
100072504	PIN, SHEAR, SLURRY CUP, HPHT CONSISTOMETER	1
100008167	SCREW SET, 6-32 X 1/4, HEX CAP	1
100010310	GASKET CUP PLUG	1