



Fann Instrument Company

Included Parts

101439783

PORTA LAB MODEL 821 (RIG LAB) LH RIG AIR

Part Number	Description	Qty
208613	ADAPTER RIG AIR	1
205852	BRUSH LARGE 2-1/2in. DIA X 12in.	1
205854	BRUSH SMALL 3/4 X 3-1/2 X 9in.	1
207999	CAP BASE w/FILTRATE TUBE f/SERIES 300 FP	1
208608	CARTRIDGES CARBON DIOXIDE CO2 10/BOX	1
208079	CELL BODY f/SERIES 300 FILTER-PRESS	1
205629	CHAIN SASH	3
205727	CLIP FUSE PIPETTE	2
208463	CLIP GRADUATE	1
205728	CLIP LARGE	2
205726	CLIP MEDIUM	3
205725	CLIP SMALL	5
205869	CYLINDER GRADUATED GLASS 10ml TC	1
205868	CYLINDER GRADUATED GLASS 25ml TC	1
205902	DISH TITRATION PLASTIC 140ml	1
206051	FILTER PAPER API 3.5in DIA (9cm) 100/BOX	2
205591	FLANGE 1/4in. f/RIG LAB	1
209387	FRAME BARE WALL MOUNT FILTER-PRESS	1
209662	FUNNEL SAND CONTENT KIT	1
206884	FUNNEL VISCOSITY MARSH PLASTIC	1
208255	GASKET API FILTER-PRESS RUBBER	6
101472471	GASKET, FLAT, 0.38 DIA.	1
205291	HOOK 'S' 1 inch	4
208604	HOSE LOW PRESSURE AIR 1 FT	1
209389	INSERT THREADED API FILTER-PRESS	1
210387	INSTRUCTION RIG LAB MODEL 821	1
205774	KNOB 5/8in. CHROME	4
206769	MUD BALANCE MODEL 14001 without CASE	1
206889	MUD CUP, MEASURING 1000ml PLASTIC - RED	1
208059	NIPPLE 1/4 X EASTMAN	2
210053	pH PAPER DISPENSER 2 ROLL pH 2-10 & 1-11	1
210060	pH PAPER DSPENSER 2 ROLL pH10-12&12.5-14	1
207444	PIN ROLL 3/32 X 3/8 STAINLESS	1
206026	PIPETTE 1 ml SEROLOGICAL	1
206029	PIPETTE 10ml SEROLOGICAL	1
209850	POTASSIUM CHROMATE INDICATOR 2oz	1
208648	PRESSURING UNIT CO2 w/TOP CAP	1
206031	ROD STIRRING 4in. PLASTIC	1
206031	ROD STIRRING 4in. PLASTIC	1
208310	SCREEN 60 MESH f/API FILTER-PRESS	2
209661	SCREEN PLASTIC SAND CONTENT	1
205809	STOPPER #2 RUBBER w/1 HOLE	1
206898	STOPWATCH DIGITAL ELECTRONIC	1
206038	THERMOMETER DIAL 50-350 DEG F 5in. STEM	1
206036	TIMER 2-HOUR	1
209910	SILVER NITRATE SOLUTION .0282N 8oz	1
209916	SILVER NITRATE SOLUTION .282N 8oz	1



**Fann Instrument Company**

**Included Parts**

208424	T-SCREW FRAME SERIES 300 FILTER-PRESS	1
209663	TUBE SAND CONTENT 100ml GLASS	1
209944	WATER DISTILLED 8oz	1