

Part No.102416675
HPHT SafeCell – Single Ended – w/CellTell
Filter Disc

Item No.	Part No.	Quantity	Description
1	102313352	1	CAP RETAINER, THREADED, MONEL K500
2	102416736	1	CELL CAP, 17-4 SS
3	102416734	1	CELL BODY, THREADED, 17-4 SS, W/ CELLTELL™
4	210538	1	FILTER DISC, CERAMIC, API-20, 10/BOX
5	209496	2	VALVE STEM, 416 SS, HARDENED RC30+
6	205649	9	O-RING, 0.176 ID x 0.0707 W, VITON
7	102365275	4	O-RING, 2-1/4 x1/8, VITON

