



Fann Instrument Company

Included Parts

232208
CLAMP ASSY, HPHT CELL 500ML

Part Number	Description	Qty
347902	BEARING, HANDLE, HPHT CELL CLAMP	2
232180	BASE, HPHT CELL CLAMP	1
232182	TOP PLATE, HPHT CELL CLAMP	1
232204	HANDLE, CLAMP, SGL LEVER, STL, 5/8-11	2
232206	ALL THREAD 5/8-11 X 15 A193B7 PL	2
100031014	NUT - HEX - 5/8-11UNC-2B - ASTM-A194 GR7	8