



Fann Instrument Company

Included Parts

**206625**  
**REGULATOR HIGH PRESSURE**

<b>Part Number</b>	<b>Description</b>	<b>Qty</b>
205651	ORING 5/16 X 1/16 NITRILE B-46 #01	2
208612	BARREL CO2 CARTRIDGE	1
206605	REGULATOR LOW PRESSURE	1
206608	FITTING TOP VALVE HTHP	1
206609	RING LOCK REGULATOR f/107C FILTER-PRESS	1
208614	ADAPTER HEAD CO2 CARTRIDGE	1
208655	VALVE NEEDLE 1/4 X 1/4 CHROME	1
209466	GAUGE 1000 PSI 2in. DIAL 1/4 BOTTOM CONN	1
209468	GAUGE 1500 PSI 2in. DIAL 1/4 BOTTOM CONN	1