



Fann Instrument Company

Included Parts

**206845**  
**PERMEABILITY PLUGGING APPARATUS 5,000 PSI 115 VOLT**

<b>Part Number</b>	<b>Description</b>	<b>Qty</b>
209471	PRESSURING UNIT CO2 HPHT FLTRP	1
209445	RECEIVER ASSEMBLY BACK PRES 100ml f/PPA	1
206850	CELL ASSEMBLY PPA 5,000 PSIG 1/4in. DISC	1
209437	HEATING JACKET f/PPA 115V	1
206869	PUMP ASSEMBLY PPA 5,000 PSIG	1
206039	THERMOMETER DIAL 50-500 DEG F 8in. STEM	1
206870	BRACKET ASSEMBLY GRADUATE CYLINDER	1
205622	TUBING TYGON 3/16 X 1/16	2
205868	CYLINDER GRADUATED GLASS 25ml TC	1
210538	FILTER DISC CERAMIC 10 MICRON 10/BOX	1
208608	CARTRIDGES CARBON DIOXIDE CO2 10/BOX	1
205869	CYLINDER GRADUATED GLASS 10ml TC	1