



Fann Instrument Company

Included Parts

206923
EP/LUBRICITY TESTER MODEL 212

Part Number	Description	Qty
205688	BEARING BALL DOUBLE SEALED	1
205687	BEARING ROLLER LOWER EP/LUB TESTER	1
206938	BELT GUARD MODEL 212 EP/LUB TESTER	1
206939	BRACKET REED SWITCH MODEL 212	1
206942	BUSHING/NUT TORQUE SET f/MODEL 212	1
205768	CABLE POWER 115V 16 AWG MALE PLU	1
205781	CIRCUIT BOARD MODEL 212 EP/LUB TESTER	1
205747	CIRCUIT BREAKER EP/LUB TESTER	1
204336	CLAMP 3/8in. DIA CABLE	1
204303	CLAMP CABLE 3/8in.	3
206932	CLAMP TORQUE WRENCH MODEL 212	1
205713	CONNECTOR ANGLE EP&LUB TESTER	1
205712	CONNECTOR MOLEX 2 PIN HSG w/.1	1
205786	CONTROLLER MOTOR SPEED f/MODEL 212	1
206896	CUP SAMPLE f/EP/LUBRICITY TESTER	1
204271	DISCONNECT FEMALE 16-14 AWG	7
204270	DISCONNECT FEMALE 90DEG FLAG 1	6
205779	FEET RUBBER 3/4in.	4
205626	GREASE ZERK 1/4-28	1
206934	HANDLE TORQUE ADJUSTMENT f/MODEL 212	1
206929	HOLDER TEST BLOCK f/MODEL 212	1
206905	INSTRUCTION EP/LUBE TESTER	1
204209	INSTRUCTION RING-BLOCK PAIRS	1
205278	KEY WOODRUFF #5 STEEL MODEL 21	1
101453815	KNOB EP/LUB TESTER (SET)	1
206944	MAGNET f/MODEL 212	2
205763	METER MODEL 212 EP/LUB TESTER	1
205785	MOTOR MODEL 212 EP/LUB TESTER	1
205678	NUT LOCK BEARING	1
206931	NUT TEST RING f/EP/LUBRICITY TESTER	1
206933	NUT TORQUE ADJUSTMENT f/EP/LUB TESTER	1
206940	PANEL MODEL 212	1
203696	PIN CONNECTOR MOLEX	2
206950	PLATE BOTTOM MODEL 212	1
205704	POTENTIOMETER RPM TORQUE MODEL 212	1
206926	RETAINER TOP BEARING MODEL 212	1
205686	RING INNER LOWER BEARING	1
205655	RING QUAD 1/2in. ID	1
205696	SEAL LOWER SHAFT MODEL 212	1
205695	SEAL UPPER SHAFT MODEL 212	1
206927	SHAFT LOCK MODEL 212	1
206930	SHAFT TEST RING MODEL 212 EP/LUB TESTER	1
206935	SHAFT TORQUE ADJUSTMENT MODEL 212	1
206928	SHAFT TORQUE MODEL 212 EP/LUB TESTER	1
205675	SHEAVE & HUB DRIVEN EP/LUB TESTER	1
205674	SHEAVE & HUB MOTOR EP/LUB TESTER	1
206941	SPACER BELT GUARD MODEL 212	2



Fann Instrument Company

Included Parts

206937	STAND CUP MODEL 212	1
206936	SUPPORT CUP STAND MODEL 212	1
205749	SWITCH REED	1
205748	SWITCH RPM-TORQUE DPDT	1
206901	TEST BLOCK EP TEST	1
206903	TEST BLOCK LUBRICITY TEST	1
206900	TEST RING EP TEST	1
206902	TEST RING LUBRICITY TEST	1
208349	THUMB-SCREW 1/4-20 X 3/4 STAINLESS	1
205670	V-BELT MODEL 212 EP/LUBRICITY TESTER	1
205677	WASHER LOCK BEARING	1
206946	WASHER TORQUE ADJUSTMENT	1
205641	WRENCH HEX KEY 3/16 SHORT ARM	1
205637	WRENCH TEST RING NUT	1
205638	WRENCH TORQUE f/EP/LUBE TESTER	1
206947	WRENCHES TORQUE BUSHING 212	1