



Fann Instrument Company

Included Parts

**206966**  
**CUP-HEATER RHEOMETER 230 VOLT**

<b>Part Number</b>	<b>Description</b>	<b>Qty</b>
206962	BODY CUP HEATER MODEL 25120	1
208865	CABLE POWER 230V 18 AWG M&F PLUG	1
209523	CAPACITOR 10000 PF 1000V DISC	1
206963	COVER CUP HEATER MODEL 25120	1
207030	CUP SAMPLE f/MODEL 280 & 286 RHEOMETER	1
205734	HEATER ROD CUP-HEATER 50 WATT 115 VOLT	2
205777	KNOB ROUND BLACK ABS w/BRASS BUSHING	1
209488	LAMP NEON HPHT FP & CUP-HEATER	1
206964	PIN-POSITIONING MODEL 25120	3
209491	PLATE THERMOSTAT MOUNT	1
205706	RESISTOR 5 OHM 1W	1
205767	STRAIN RELIEF .325-.360 DIA-.0	1
205798	THERMOSTAT 50-300 DEG F HEATING CUP	1