



Fann Instrument Company

Included Parts

207083
CALIBRATION STAND f/280 & 286 RHEOMETER

Part Number	Description	Qty
205052	WEIGHT 20gm BRASS WITH HOOK	1
205053	WEIGHT 50gm BRASS WITH HOOK	1
205683	BEARING BOB-SHAFT 280 & 286 RHEOMETER	1
205778	FEET RUBBER 1/2in.	4
206987	JEWEL MODEL 280 286 RHEOMETER	1
207849	THREAD RED CALIBRATION KITS	1