



Fann Instrument Company

Included Parts

**207205**  
**GEAR BOX MODEL SR12 50Hz**

<b>Part Number</b>	<b>Description</b>	<b>Qty</b>
207480	GEAR BOX SR-12 & SASR-12	1
207203	PLATE CHANGE GEAR 35SA/SR-12 VISCOMETER	1
207475	POST #1 GEAR SR-12 & SA/SR-12	1
207476	POST #2 GEAR SR-12 & SA/SR-12	1
207814	RETAINER EXTERNAL VISCOMETER & 12BL	2
207478	GEAR IDLER #1 29T-48P SR12	1
207577	GEAR ASSEMBLY IDLER #2 SA-SR12 VISCOMETR	1
207579	GEAR ASSEMBLY OUTPUT 35SA/SR12	1
207725	BEARING SR-12 GEARBOX AND VISCOMETER	2
207860	RETAINER EXTERNAL 35SA/SR12	1
207477	ARM GEAR SHIFT MODEL 35/SR12 VISCOMETER	1
207510	KNOB RED f/MODEL 35 VISCOMETER	1
207485	WASHER FLAT STEEL 1/4X1/2X1/16	2
207843	1/4-28 X 2-1/2 SHCS STAINLESS	4
208847	TAG CAUTION GEAR SHIFT SR12 VI	1