



Fann Instrument Company

Included Parts

207673
FILTER-PRESS API 6-UNIT w/MANIFOLD

Part Number	Description	Qty
207999	CAP BASE w/FILTRATE TUBE f/SERIES 300 FP	6
208062	CAP TOP SERIES 300 (API) FILTER-PRESS	6
208079	CELL BODY f/SERIES 300 FILTER-PRESS	6
205587	COUPLING HEX 1/4 FNPT STAINLESS	6
205869	CYLINDER GRADUATED GLASS 10ml TC	6
206051	FILTER PAPER API 3.5in DIA (9cm) 100/BOX	3
207674	FRAME 6-UNIT FILTER-PRESS	1
208255	GASKET API FILTER-PRESS RUBBER	36
208604	HOSE LOW PRESSURE AIR 1 FT	6
209389	INSERT THREADED API FILTER-PRESS	6
207128	INSTRUCTION FILTER-PRESS SERIES 300	1
207675	MANIFOLD 6-UNIT FILTER-PRESS	1
208059	NIPPLE 1/4 X EASTMAN	7
207444	PIN ROLL 3/32 X 3/8 STAINLESS	6
205589	PLUG PIPE 1/4 NPT BRASS	1
208283	ROD SUPPORT API FILTER-PRESS	6
208310	SCREEN 60 MESH f/API FILTER-PRESS	6
208626	STOP MUD FELT 33700 REGULATOR	6
208459	SUPPORT GRADUATE	6
208349	THUMB-SCREW 1/4-20 X 3/4 STAINLESS	6
207678	TRAY CELL PREPARATION 6 UNIT	1
208424	T-SCREW FRAME SERIES 300 FILTER-PRESS	6
208653	VALVE BLEEDER 1/4in. NPT CHROME PLATED	6
208655	VALVE NEEDLE 1/4 X 1/4 CHROME	6
204218	WARRANTY CARD	1