



Fann Instrument Company

Included Parts

208594
DEAD WEIGHT HYDRAULIC ASSEMBLY

Part Number	Description	Qty
204211	INSTRUCTION DEADWEIGHT HYDRAULIC FP	1
205653	ORING 1/2 X 3/32 NITRILE B-46 #112	1
207889	WRENCH HEX KEY 1/8 SHORT ARM	1
207929	HOSE LOW PRESSURE AIR 3 FT	1
208599	VALVE CHECK 1/8 MNPT STAINLESS STEEL	2
208602	GAUGE 160 PSI 1.5in DIAL 1/8 BOTTOM CONN	1
208603	ADAPTER 1/4 X 1/8 NPT SPECIAL TAPER	1
208759	VALVE NEEDLE 1/8 X 1/8 CHROME	1