



Fann Instrument Company

Included Parts

208760
MIXER LAB MODEL 230V w/STAND

| Part Number | Description | Qty |
|--------------------|--------------------------------------|------------|
| 206537 | BASE STAND f/LAB MIXER | 1 |
| 206560 | IMPELLER LAB & PORTABLE MIXERS | 1 |
| 206595 | ARM MOTOR LABORATORY MIXER | 1 |
| 206596 | KNOB ASSEMBLY 5/16-18 AL LAB MXR | 1 |
| 206597 | DAMPER SLIDE LABORATORY MIXER | 1 |
| 206664 | CABLE POWER RETRACTABLE MIXER | 1 |
| 208761 | MOTOR 230V LAB MIXER 35301 | 1 |
| 206733 | SCREW MOTOR ARM TO BRACKET LAB MIXER | 3 |
| 206594 | BOX JUNCTION f/LAB MIXER | 1 |
| 205767 | STRAIN RELIEF .325-.360 DIA-.0 | 1 |
| 205718 | PLUG 2 POLE 3 WIRE 15 AMP 250 | 1 |