



Fann Instrument Company

Included Parts

209529

FILTER-PRESS HPHT 175ml 230V DBL END N2O

Part Number	Description	Qty
209502	RECEIVER ASSEMBLY BACK PRESS 15ml N2	1
209506	HEATING JACKET f/175ml 230 VOLT	1
209545	MANIFOLD DUAL HPHT FILTERPRESS	1
209486	INSTRUCTION FILTER-PRESS HPHT SERIES 387	1
208654	WRENCH ADJUSTABLE 6in.	1
206056	FILTER PAPER 2.5in. DIA (6.35cm) 100/BOX	1
205649	ORING 3/16 X 1/16 VITON #008-7	12
205662	ORING 2-1/4 X 1/8 NITRILE B46 #228	3
206038	THERMOMETER DIAL 50-350 DEG F 5in. STEM	1
205868	CYLINDER GRADUATED GLASS 25ml TC	1
209530	CELL HPHT 175ml 303SS 1800PSI DOUBLE END	1