



Fann Instrument Company

Included Parts

209531
CELL HPHT 175ml 303SS 1800PSI SINGLE END

Part Number	Description	Qty
209501	SCREW CAP LOCKING 5/16-18	8
209536	CAP HPHT 1800 PSIG DETACHABLE SCREEN	1
209496	VALVE-STEM HPHT FILTER-PRESS	2
205649	ORING 3/16 X 1/16 VITON #008-7	8
205662	ORING 2-1/4 X 1/8 NITRILE B46 #228	2
209534	SCREEN 325 w/60 MESH BACK-UP	1
209493	WRENCH HEX KEY 5/32	1
209569	CELL BODY 175ml 303SS w/PRESS INDICATOR	1