



Fann Instrument Company

Included Parts

209545
MANIFOLD DUAL HPHT FILTERPRESS

Part Number	Description	Qty
208759	VALVE NEEDLE 1/8 X 1/8 CHROME	1
209466	GAUGE 1000 PSI 2in. DIAL 1/4 BOTTOM CONN	1
209467	GAUGE 3000 PSI 2in. DIAL 1/4 BOTTOM CONN	1
209470	REGULATOR 1500 PSI	1
209473	GLAND LONG NITROGEN SERVICE	1
209474	HOSE 3000 PSI 3 FT X 3/16in.	2
209475	BODY DUAL MANIFOLD f/MODEL 38765	1
209476	VALVE RELIEF 1/8in. 1200 PSI	1
209477	VALVE RELIEF 1/8in. 750 PSI	1
209515	BLOCK MANIFOLD HPHT FILTER-PRESS	1
204648	PIN COX SAFETY LOCKING	1
204894	REGULATOR MODIFIED 3841141003	1
205583	NIPPLE 1/4 NPT HEX STAINLESS	2
205585	TEE STREET 1/8in. STAINLESS STEEL	2
205589	PLUG PIPE 1/4 NPT BRASS	1
205595	UNION HALF 1/4 X 1/8 PLATED	1
205596	CONNECTOR MALE 1/4 NPT X 45 STAINLESS	1
205597	ELL 1/8FNPT X 1/4 37D FLARE STAINLESS	1
203950	NUT GLAND RIGHT HAND	1
209472	NUT GLAND LEFT HAND	1