



Fann Instrument Company

Included Parts

**209558**

***FILTER-PRESS HPHT 500ml 230V DBL END N2O***

<b>Part Number</b>	<b>Description</b>	<b>Qty</b>
209553	BASE MODEL 38770 HPHT FILTER-PRESS	1
208865	CABLE POWER 230V 18 AWG M&F PLUG	1
209523	CAPACITOR 10000 PF 1000V DISC	1
209552	CELL HPHT 500ml 303SS 1800PSI DOUBLE END	1
209556	COVER HPHT FILTER-PRESS 38770	1
205868	CYLINDER GRADUATED GLASS 25ml TC	1
205779	FEET RUBBER 3/4in.	4
206058	FILTER DYNALLOY X-5 2.5in DIA (6.35 CM)	1
206056	FILTER PAPER 2.5in. DIA (6.35cm) 100/BOX	1
205736	HEATER ROD HPHT JACKET 200 WATT 115 VOLT	4
209486	INSTRUCTION FILTER-PRESS HPHT SERIES 387	1
205393	INSULATION HEATER 1/32 THK J50	1
205777	KNOB ROUND BLACK ABS w/BRASS BUSHING	1
209488	LAMP NEON HPHT	1
209554	LEG HPHT FILTER-PRESS SERIES 38770	2
209545	MANIFOLD DUAL HPHT FILTERPRESS	1
205668	ORING 1-3/8 X 1/16 NITRILE B-46 #0	6
205662	ORING 2-1/4 X 1/8 NITRILE B46 #228	24
205649	ORING 3/16 X 1/16 VITON #008-7	24
209491	PLATE THERMOSTAT MOUNT	1
209542	RECEIVER ASSEMBLY BACK PRESS 100ml N2	1
205706	RESISTOR 5 OHM 1W	1
209432	SCREW PIN LOCATOR HEATING WELL	1
206039	THERMOMETER DIAL 50-500 DEG F 8in. STEM	1
205801	THERMOSTAT 50-600 DEG F f/HPHT FP	1
205595	UNION HALF 1/4 X 1/8 PLATED	1
209555	WELL HPHT FILTER PRESS 38770	1
208654	WRENCH ADJUSTABLE 6in.	1