

209694
KIT METHYLENE BLUE TEST 230V

Part Number	Description	Qty
203397	10-32 X 3/4 RHMS STAINLESS	2
203421	10-32 X 1/4 RHMS STAINLESS	3
203422	8 X 3/4 PHSMS STAINLESS	2
205283	BAG PLASTIC ZIPLOCK 4in. X 6in.	1
205290	STOPPER #6 RUBBER	1
205726	CLIP MEDIUM	1
205728	CLIP LARGE	1
205729	CLIP FLASK	1
205742	HOT PLATE 230V 300 WATT w/THERMOSTAT	1
205774	KNOB 5/8in. CHROME	1
205799	NEEDLE LUER-LOK SYRINGE	1
205847	BOTTLE 2 oz AMBER DROPPER GLASS 20-400	1
205867	CYLINDER GRADUATED GLASS 50ml TD	1
205895	SYRINGE 5cc LUER-LOCK GLASS	1
205914	FLASK ERLNMEYER GLASS 250ml	1
206026	PIPETTE 1 ml SEROLOGICAL	1
206029	PIPETTE 10ml SEROLOGICAL	1
206030	ROD STIRRING 10in. GLASS	1
206730	PUMP PIPETTE FILLER/DISPENSER	1
207633	NUT 10-32 HEX REGULAR STAINLESS	2
207759	8-32 X 1/4 BHMS STAINLESS	1
208704	WASHER SPLIT #10 STAINLESS STEEL	2
209680	INSTRUCTION METHYLENE BLUE KIT	1
209682	TAG NAMEPLATE METHYLENE BLUE	1
209687	METHYLENE BLUE .01 meq 8oz	1
209796	BRACKET f/HANDLE No. 50003	1
209797	HANDLE PLASTIC BLACK	1
209798	CASE MULTI-KIT STAINLESS STEEL	1
209944	WATER DISTILLED 8oz	1
210070	PIPETTE DROPPER 1 ml	1
210071	BULB RUBBER SUCTION f/653 SAMPLE CELL	1
209847	HYDROGEN PEROXIDE 3% 8oz	1
206050	FILTER PAPER 5in. DIA (12.5cm) 100/BOX	1

Item not included in Kit - Required for API Test - Order Separately

209874 SULFURIC ACID 5N 8oz