



Fann Instrument Company

Included Parts

**210288**

**CELL AGING NON-PRESSURIZED 260 ml 303 STAINLESS**

<b>Part Number</b>	<b>Description</b>	<b>Qty</b>
210307	CAP INNER AGING CELL w/o VALVE	1
210306	CAP OUTER AGING CELL w/o VALVE	1
210308	CELL BODY AGING CELL 260 ml 303 STAINLESS	1
210276	GASKET AGING CELL TEFLON	1
210287	INSTRUCTION AGING CELL	1
210278	SCREW SET 5/8 in	1
210277	WRENCH HEX KEY 5/16	1