



Fann Instrument Company

Included Parts

210294
CELL CORROSION TEST 500ml 303SS

Part Number	Description	Qty
210293	CAP INNER CORROSION TEST CELL	1
210311	CAP OUTER AGING CELL	1
210309	CELL BODY 500ml 303 STAINLESS	1
210276	GASKET AGING CELL TEFLON	1
210086	GROMMETS & SCREWS COUPON HOLDER	1
210287	INSTRUCTION AGING CELL	1
205649	ORING 3/16 X 1/16 VITON #008-7	2
210281	SCREW SET 3/8in.	3
210282	VALVE-STEM 416 STAINLESS	1
210284	WASHER FLAT GALVANIZED	1
205641	WRENCH HEX KEY 3/16	1