

**210465**  
**KIT RETORT 50ml OIL & WATER 115 VAC**

<b>Part Number</b>	<b>Description</b>	<b>Qty</b>
205855	BRUSH TEST TUBE	1
205850	BRUSH WIRE-ALL STAINLESS STEEL 1/8 DIA	1
205769	CABLE EXTENSION 10' 115VAC	1
210467	CASE 50ml RETORT STAINLESS STEEL	1
210441	CLEANER PIPE 5mm 20/PKG	1
205728	CLIP LARGE	3
205726	CLIP MEDIUM	2
204295	CLOSURE SNAP 5/8	1
210439	CORKSCREW	1
208776	CYLINDER GRADUATED GLASS 50ml TC	1
204271	DISCONNECT FEMALE 16-14 AWG	2
206118	DRILL TWIST (#31) 6 IN. LENGTH	1
210474	HEATER BLOCK RETORT 50ml 115/230 V	1
205737	HEATER ROD 871 RETORT 350 WATT 120 VOLT	2
205758	HOLDER LAMP OVEN & RETORT 115 VOLT	1
210466	INSTRUCTION RETORT KIT 50ml	1
210459	INSULATION BULK	20
205759	LAMP NEON OVEN & RETORT 115 VOLT	1
205757	LENS CLEAR 10ML RETORT 12 VOLT	1
210435	LUBRICANT HIGH TEMPERATURE 1oz TUBE	1
210468	MACHINE PARTS SET 50ml RETORT	1
205744	PLATE SWITCH ON / OFF	1
205719	RECEPTACLE RECESSED MALE 2 PRO	1
210446	SPACER RECEPTACLE	1
210433	SPATULA RETORT	1
210440	STEEL WOOL 1/4 POUND PACKAGE	1

**fann**<sup>®</sup>  
Fann Instrument Company  
Included Parts

**fann**<sup>®</sup>

**Fann Instrument Company**

**Included Parts**

205147

TERMINAL STRIP CINCH JONES 2-1

1

**fann**<sup>®</sup>  
Fann Instrument Company  
Included Parts

**fann**<sup>®</sup>  
Fann Instrument Company  
Included Parts

**fann**<sup>®</sup>

**Fann Instrument Company**

**Included Parts**

206119

WRENCH TAP T-HANDLE

1