



Fann Instrument Company

Included Parts

210528

MULTI-RETORT 5 CHAMBER 230 VOLT

Part Number	Description	Qty
208627	BALL 1/4in. DIAMETER STAINLESS STEEL	1
205853	BRUSH MEDIUM 1-1/2in. X 10-3/4in.	1
205770	CABLE 14 AWG 3 CONDUCTOR SJO 300	8
204303	CLAMP CABLE 3/8in.	4
210441	CLEANER PIPE 5mm 20/PKG	1
205728	CLIP LARGE	5
204339	CONTROLLER PROPORTIONAL SOLID	1
210439	CORKSCREW	1
205867	CYLINDER GRADUATED GLASS 50ml TD	5
204273	DISCONNECT 2 MALE TABS	10
204344	DISCONNECT FEMALE .250 TAB 22-	2
204271	DISCONNECT FEMALE 16-14 AWG	48
204270	DISCONNECT FEMALE 90DEG FLAG 1	2
339350	FAN 3.125in. SQ X 1.5in. THK 2305V 37CFM	2
205779	FEET RUBBER 3/4in.	6
210534	FUSE LINK-TEMP MULTI-RETORT	1
205391	HEATER ROD MODEL 90 & MULTI-RETORT	10
210533	INSTRUCTION MULTI-RETORT	1
210530	LIFTER TOOL MULTI-RETORT	1
210435	LUBRICANT HIGH TEMPERATURE	1
205648	ORING 5/32 X 1/16 NITRILE B-46 #00	5
205718	PLUG 2 POLE 3 WIRE 15 AMP 250	1
204667	POTENTIOMETER 250K OHM	1
204903	RELAY 240 VOLT 30 AMP DPDT	1
204927	RETAINER INTERNAL 1in. BORES	1
210535	RETORT CHAMBER 50ml MULTI-RETORT	5
203445	SPACER 1/4 OD X #6 CLEAR ID X	4
210433	SPATULA RETORT	1
207009	SPRING	1
210440	STEEL WOOL 1/4 POUND PACKAGE	2
210529	SWITCH OVER-TEMP MULTI RETORT	1
205122	SWITCH ROCKER BREAKER 15 A	1
208801	THERMOCOUPLE TJ72-ICSS-18U-6304 SS	1
208030	TUBE FILTRATE DRAIN	5
209938	WETTING AGENT 1oz	1