



Fann Instrument Company

Included Parts

381464

CHILLER ASSY MOD 75 VISCOMETER 220/50 HZ

Part Number	Description	Qty
204082	ELBOW 1/4" MNPT X 1/4" TUBE BRASS	2
205584	TEE 1/4" BRASS	2
208651	GAUGE 200 PSI 2" DIAL, 1/4" BOTTOM CONNECTION	1
204085	CONNECTOR MALE 1/4" TUBE X 1/4" MNPT	1
204083	TEE MALE BRANCH 1/4" MNPT X 1/4" TUBE	2
205405	TUBING, NYLON FLEXIBLE, 1/4" OD	12
206158	VALVE NEEDLE 1/4" FNPT VITON	1
204076	INSULATION RUBBER TUBE 3/8" ID	18
204077	INSULATION RUBBER TUBE 3/4" ID	8
204075	COUPLING QUICK-INSERT 1/4" FNPT	1
204074	COUPLING QUICK-DISCONNECT 1/4" MNPT BRASS	1
381403	CHILLER RECIRCULATING 220 VOLT 50Hz	1
204875	ETHYLENE GLYCOL INDICATOR GRADE 32oz	2
209947	WATER DISTILLED 1 GALLON	1
204080	CAP CELL INSULATOR	1
205589	PLUG PIPE 1/4" NPT BRASS	1
208977	TUBING 1/4" OD X .035 WALL 304 SS	8
208854	CAP PIPE 1/4" NPT	1