

101764749
SPARES 1 YEAR RHEOMETER MODEL iX77

Part Number	Description	Qty
100123611	SCREW, SET, #6-32 X 3/16, ALLEN, KCP, SS1 (UNBRAKO)	20
101623710	O-RING, SIZE 2-226, COLOR WHITE, 600 DEG F, ID 1.984 IN X .139 IN WALL	6
101767884	O-RING 2-153 VITON 75 DURO 3.487 ID	6
101626358	FILTER, 40 MICRON, 304 WIRE MESH, 3B FILTER	4
101626628	BUSHING, ROTOR, BRONZE	3
101629340	KIT, SOFT GOODS, 30,000 PSIG, TESCOM REGULATOR	3
101629342	KIT, REPAIR, TESCOM 54-21N0T44A-074 REGULATOR, 30,000 PSIG	1
101630287	KIT, 3 O-RINGS f/P8791-HP-40WEV ORING 212V, 138V, AND 138VP	3
101736752	O-RING, SIZE 2-220, COLOR WHITE, 600 DEG F, .139 IN WALL	6
203365	VEE JEWEL	5
203366	BALL 3/16in. DIAMETER STAINLESS STEEL	5
203367	BEARING ANGULAR CONTACT MODEL iX77	3
203509	BOTTLE 16oz CLEAR PET 63mm	2
204607	O-RING 0.209 ID X .070 NITRILE C67-90	10
204816	GREASE HIGH TEMP PERMATEX 82325	2
207122	CALIBRATION FLUID 200 cP 16oz	4
208939	RING BACKUP TOP CELL	5
208940	RING BACKUP BOTTOM CELL	5
210435	LUBRICANT HIGH TEMPERATURE 1 OUNCE TUBE	3