



Fann Instrument Company

Spare Parts

101877938

SPARES 1 YEAR MACS II CEMENT ANALYZER

Part Number	Description	Qty
101667558	THERMOCOUPLE, TOP	1
101671856	THERMOCOUPLE, BOTTOM	1
101918773	THERMOCOUPLE, SAFETY	1
101594774	BACKUP RING, PRESSURE VESSSEL	8
101710014	O-RING, ULTRA PAROFLUOR, ID 2.984, WHITE	4
101710015	O-RING, ULTRA PAROFLUOR, ID .739, WHITE	4
101596797	TIMING BELT 3/8W X 120 T	2
101542809	RUPTURE DISC, BURST PRESS 34500 PSI @ 72 F	1
101636441	PADDLE ASSEMBLY	5
100001522	RING RETAINING EXTERNAL 0.25	50
101836939	GASKET, HIGH TEMP, CAP PLUG	30
101836938	GASKET, HIGH TEMP, CUP BASE	10
101837384	O-RING, 568-232, PERFLUORO WHITE 90 DUR 600F	4
101626358	FILTER ELEMENT, 25 MICRON, HP 304 SS	1
101629340	KIT, SOFT GOODS, 30 KSI TESCOM REGULATOR	3
101629342	KIT, REPAIR, TESCOM 54-21N0T44A-095	3
101837396	O-RING, 568-107, PERFLUORO WHITE 90 DUR 600F	4
100002427	O-RING, TEFLON, 3/8OD X 1/4ID	50
102052994	WATER FILTER ELEMENT, 10 MICRON	5
101673153	TUBING, 1/4 IN, POLYURETHANE	120"
101485964	TUBING, 3/8 IN O.D. POLYURETHANE	10'
101630287	KIT,3 ORINGS FOR P8791-HP-40WEV ORING 2	4