



Fann Instrument Company

Spare Parts

**206333**  
**SPARES 1 YEAR STICKING TESTER**

<b>Part Number</b>	<b>Description</b>	<b>Qty</b>
204648	PIN COX SAFETY LOCKING	1
205649	ORING 3/16 X 1/16 VITON #008-7	20
205652	ORING 3/8 X 1/16 NITRILE B-46 #012	8
205662	ORING 2-1/4 X 1/8 NITRILE B46 #228	2
205665	ORING 3-3/4 X 1/8 NITRILE B-46 #23	2
206056	FILTER PAPER 2.5in. DIA (6.35cm) 100/BOX	1
206908	GASKET STICKING TESTER NEOPRENE	10
206909	GASKET STICKING TESTER PLASTIC	4
206910	DISC LOCKING MESH f/STICKING TESTER	10
207232	SCREEN 60 MESH f/HPHT & 1/2 AREA FLTRPRS	4
208608	CARTRIDGES CARBON DIOXIDE CO2 10/BOX	20
208616	SEAT CO2 REGULATOR 33700	1
208626	STOP MUD FELT 33700 REGULATOR	6
208643	KIT REGULATOR REPAIR	1
209496	VALVE-STEM HPHT FILTER-PRESS	2
101472471	GASKET, FLAT, 0.38 DIA.	12