



Fann Instrument Company

Spare Parts

**206360**  
**SPARES 1 YEAR VISCOMETER MODEL 70/75**

<b>Part Number</b>	<b>Description</b>	<b>Qty</b>
203365	VEE JWLC, CELL ASSY, M70, 75 & 77 VISMR	2
203367	BEARING ANGULAR CONTACT MODEL 70 & 75	24
203994	COUPLING MALE 3/8NPT X 1/4TUBE	1
203996	STEM VALVE VENT STAINLESS STEEL M70	2
203998	PACKING VALVE VENT TEFLON M70	2
204604	ORING 1.359 X .139 NITRILE C67-90	60
204605	ORING 2 X 1/8 VITON V14 #226-75	120
204607	ORING 0.209 ID X .070 NITRILE C67-90	36
205203	KIT PUMP REPAIR MODEL 70 (ENERPAC)	1
206092	TRANSISTOR POWER MOSFET MTP12N	1
207708	DISC RUPTURE 23,000 PSI 1/4 INCONEL	1
207874	OIL HEAD MODEL 50 RHEOMETER 1 QT	12
208898	HINGE PIVOT BUSHING M70/75 VISCOMETER	6
208905	BSHG, RTR M70, 75 & 77 VISMR	12
208939	RING BACKUP TOP CELL MODEL 70	1
208940	RING BACKUP BOTTOM CELL MODEL 70 & 75	1