



Fann Instrument Company

Spare Parts

206531
SPARES 2 YEAR RETORT KIT 230 VOLT

Part Number	Description	Qty
205716	CONNECTOR MALE RECESSED WITH G	1
205738	HEATER ROD 871 RETORT 350 WATT 230 VOLT	1
205755	SWITCH THERMOSTAT L8262	1
205758	HOLDER LAMP OVEN & RETORT 115 VOLT	2
205759	LAMP NEON OVEN & RETORT 115 VOLT	2
205766	CABLE ASSEMBLY 230 VOLT POWER	1
205802	THERMOSTAT 1100 DEG F f/RETORTS	1
205854	BRUSH SMALL 3/4 X 3-1/2 X 9in.	4
205869	CYLINDER GRADUATED GLASS 10ml TC	8
209938	WETTING AGENT 1oz	6
210433	SPATULA RETORT	1
210434	RETORT CHAMBER 10ml	1
210435	LUBRICANT HIGH TEMPERATURE	8
210436	LID CHAMBER 10ml RETORT	4
210437	CONDENSER 10ml RETORT	1
210439	CORKSCREW	1
210440	STEEL WOOL 1/4 POUND PACKAGE	8
210441	CLEANER PIPE 5mm 20/PKG	8
210445	INSTRUCTION RETORT KIT 10ml	1